



CHASE RESEARCH CRYOGENICS LTD.
WORLD LEADERS IN SUB-KELVIN CRYOGENICS

TWO-STAGE SUB-KELVIN ^3He CRYOCOOLER

TYPE GL7 (Helium 7)

GENERIC INSTALLATION AND OPERATING INSTRUCTIONS



Photo shows a typical CRC GL7 (Helium 7) cryocooler

CRC Ltd, Cool Works, Unit 2 Neepsend Ind. Est., Parkwood Rd, Sheffield S3 8AG, UK.

Tel: +44114478 4140

Director: Dr. S.T. Chase

Registered in England & Wales, No. 4643351

Secretary: Dr. L.C. Kenny

VAT registration No. GB 763 8558 84

Contents

1. GENERAL HANDLING.....	3
1. SAFETY OF CHASE RESEARCH CRYOGENICS PRODUCTS.....	3
1.1. Pressure Equipment Directive 97/23/EC (Pressure Equipment Regulations 1999).....	3
1.2. Pressure Systems Safety Regulations 2000.....	3
1.3. Safe Operation.....	4
1.4. Risk Assessment	4
3. A BRIEF DESCRIPTION OF THE CRYOCOOLER UNIT	5
4. INSTALLATION.....	6
4.1. Mechanical	6
4.2. Electrical	6
5. ATTACHING YOUR EXPERIMENT TO THE CRYOCOOLER.....	7
5.1. Radiation shielding	8
6. OPERATION.....	8
6.1. Pre-cool	8
6.2. Running the cryocooler	10
6.3. Typical re-cycle from cold	11
7. OPTIMISING THE PERFORMANCE OF YOUR CRYOCOOLER	12
7.1. Pre-cool to 4K.....	12
7.1. Running the cryocooler	13
8. STANDARD Pin-Out Assignments.	15

THIS GENERIC OPERATING MANUAL describes how to install and operate a CRC GL7-Helium 7 cryocooler. It is accompanied by an Excel file that contains the validation test data and the calibration files that are **specific** to the cryocooler unit that you have purchased.

You are advised to make a note below of the location of the Excel file specific to your cryocooler unit.

1. GENERAL HANDLING

WARNING!

CRC CRYOCOOLERS CONTAIN HELIUM GAS AT HIGH PRESSURE.

Do not crush, twist or bend the unit. Avoid applying mechanical stresses. Do not heat the unit above room temperature. Keep in a sealed cryostat, or in the shipping box and brace in which it came.

Do not hold or lift the unit by means of the cold heads.

Do not tamper with the copper capillary fill tubes.

Avoid the use of acid fluxes when soldering in the vicinity of the cryocooler. Chloride based fluxes will corrode stainless steel and could damage your cryocooler.

After unpacking the cryocooler according to the instructions supplied, the cryocooler should be immediately transferred into the host cryostat. The shipping brace doubles as a stand for the cryocooler, though when used as a stand, the three screws through the aluminium plate into the cold heads should NOT be in place. When picking the cryocooler up, it should be firmly grasped by the cryopump radiation shield or the main plate/angle bracket.

1. SAFETY OF CHASE RESEARCH CRYOGENICS PRODUCTS

1.1. Pressure Equipment Directive 97/23/EC (Pressure Equipment Regulations 1999)

This CRC cryocooler unit is manufactured in accordance with Sound Engineering Practice. The volume and gas pressure within the cryocooler are such that the equipment falls below the lower classification limit in Annex II of the Pressure Equipment Directive. Hence the requirements for Conformity Assessment do not apply and no Declaration of Conformity can be made, or CE marking applied.

The cryocooler is covered by Article 3 Paragraph 3 of the Pressure Equipment Directive, which states: "Pressure equipment and/or assemblies below or equal to the limits in sections 1.1, 1.2 and 1.3 and section 2 respectively must be designed and manufactured in accordance with the sound engineering practice of a Member State in order to ensure safe use. Pressure equipment and/or assemblies must be accompanied by adequate instructions for use and must bear markings to permit identification of the manufacturer or of his authorized representative established within the Community. Such equipment and/or assemblies must not bear the CE marking referred to in Article 15."

1.2. Pressure Systems Safety Regulations 2000

This cryocooler unit does not contain a pressure x volume product exceeding 250 bar-litres hence PSSR regulations 5(4), 8-10 and 14 do not apply. This means that the system does not require a written scheme of examination. The cryocooler is not 'mobile' in the sense intended in the PSSR hence *the owner* has duties under these regulations to ensure that a) the safe operating limits are not exceeded; b) the unit is operated in accordance with these instructions; c) the unit is returned to Chase Research Cryogenics Ltd in the event that any maintenance is required. The cryocooler contains no user-serviceable parts.

1.3. Safe Operation

The safe operating temperature range of this cryocooler is 0 to 320 K.

1.4. Risk Assessment

CRC cryocoolers contain Helium gas under pressure. The stored energy of the system is less than 50 bar litres. All system components are integrity tested during manufacture; the slightest leak will make the cryocooler lose its stored gas and cease to function. A unit that has leaked presents no risks whatever to the user; the following risk assessment applies therefore only to functional units.

Hazards and consequences

Accidental damage to the cryocooler unit could result in the sudden release of pressurised gases, causing mechanical failure of the unit and potential injury (or damage to surrounding instruments) from ejected debris.

Possible events leading to failure are: overheating of the unit, for example in a fire; dropping or crushing of the unit; twisting or bending of the gas tubes. Mechanical damage to the unit is most likely to occur during assembly of the instrument of which the cryocooler forms part.

Risks without controls in place

It is extremely unlikely that the above events will lead to danger. Chase Research Cryogenics Ltd has produced more than one hundred cryocooler units of various designs, which are in use for a range of applications worldwide. To date there has never been a sudden failure of a cryocooler unit – indicating that with normal use (including inevitable handling mishaps) the units have an excellent safety record. User experience to date shows that accidental mechanical damage to cryocooler units is likely to result in slow leaks, not sudden failures.

Controls in place

The controls that are in place to eliminate (as far as reasonably practicable) the risks arising from mechanical damage to a cryocooler unit are:

- This written instruction manual, containing warnings about the potential risks arising from damage to the unit and alerting the user to more risky operations;
- Instructions that the unit should not be used if it has been subjected to overheating, dropping, crushing, bending or twisting;
- A warning label on the transit box that the instructions should be read prior to handling the unit.


The applications for which cryocooler units are intended make it impossible to place warning labels on the unit itself. However if the cryocooler is incorporated into another instrument, that instrument should carry a warning label to alert the user that the cryocooler contains no user-serviceable parts and should not be disassembled.

Risks with controls in place

Providing users read and follow this instruction manual the risks are negligible.

3. A BRIEF DESCRIPTION OF THE CRYOCOOLER UNIT

This cryocooler unit has two cold heads on one side of the circular main plate, as can be clearly seen in the illustration below. The cryocooler consists of two modules; the ^3He module and the ^4He buffer module. Brackets connect the ^3He module gas pipe to the ^4He module head and film burner. In operation the unit will be inverted with respect to this picture, with the cold heads at the bottom.



This model of cryocooler provides three points at which heat may be extracted from a user's experiment mounted on a separate cold table. They are the ^3He cold head, the ^4He buffer head, and the ^4He film burner. There are holes tapped on each of these surfaces for thermal connections between your experiment and the cryocooler.

Each of the two cold heads is provided with a calibrated RuO_2 or Cernox thermometer sensor to monitor the temperature. These sensors are inserted into sockets machined directly into the heads. Wiring for the thermometer sensors is carried from two sets of two isolated standoffs, one set on each head.

The two cryopumps and gas-gap heat switches are on the other side of the main plate. Each cryopump has a heater element that controls the cooling cycle. Both cryopumps are provided with standard active gas-gap heat switches that are activated by $10\text{ k}\Omega$ heater resistors, and each heat switch also carries a diode thermometer. Heat straps are fitted between the heat switches and the cryopumps. All electrical connections are brought out to a connector mounted onto the main plate. Pin-outs are listed at the end of this manual.

4. INSTALLATION

4.1. Mechanical

Before installing the unit in your cryostat, be sure to remove all of the pieces of foam board packing from around the pumps, as mentioned in the unpacking instructions.

There should be no need to touch the heat switches or heat straps during installation or normal operation of the cryocooler. The heat switches can be easily damaged, and if bent or twisted are likely to fail.

This unit is designed to work in either 'wet' cryostat using liquid ^4He , or in a 'dry' cryostat, i.e. from a mechanical cryocooler head such as a pulse tube. Mounting holes are provided on the main plate for attaching the cryocooler to your cryostat cold plate. There are twelve 4.1mm diameter (M4 or UNC6 clearance) holes symmetrically distributed upon a 115 mm pitch circle around the periphery of the circular main plate. In addition to these, there is also a row of 4 x M4 or UNC6 clearance holes at $\frac{1}{2}$ " (12.7mm) centres, close to one edge of the main plate.

Always use spring washers, or suitable low expansion washers (e.g. Invar or Tungsten), under every bolt head. These will take out differential thermal contraction that might otherwise cause loosening of the bolts, and thus compromise thermal contact. Because the cooling down of the heads depends upon gas convection, the cryocooler *must* be kept close to vertical during the cooldown process.

4.2. Electrical

All electrical connections are brought out to a 25-pin MDM-SSP connector mounted onto the main plate. Pin-outs are listed at the end of this instruction manual. Voltage / current requirements for driving the heater and thermometers are summarised in the table below.

ITEM	NUMBER	IMPEDANCE/ JUNCTION VOLTAGE	VOLTAGE/ CURRENT
3-pump heater	1 off	400 Ω approx.	0 to 20V
4-pump heater	1 off	200 Ω approx.	0 to 24V
Heat switch heaters	2 off	10k Ω	0 to 5V
Diode thermometers	5 off	0.5 to 1.8V	10 μA DC
^4He -head RuO $_2$ thermometer	1 off	1k Ω to 3k Ω	1 μA max.
^3He -head RuO $_2$ thermometer	1 off	1k Ω to 7k Ω	100nA max.

Generic (i.e. standard calibration) RuO₂ sensors from Lakeshore Cryotronics are the default option on all CRC cryocoolers. Individually calibrated 'CERNOX' or RuO₂ sensors are only fitted (at additional cost) at the customer's express requirement. Generic diode calibration curves for the cryopump and heat switch diodes are supplied as standard by CRC Ltd. A calibration curve specific to the film burner diode is also supplied in the Excel data file that accompanies each cryocooler.

Wiring for the temperature sensors on the cold heads is carried to insulated stand-offs mounted on each cold head, and from there directly to the 25-pin connector. The thermometer on the ³He cold head is operated as a 4-wire device and should be excited with an AC current no greater than 100nA, corresponding to a voltage of around 2mV at base temperature. The thermometer on the ⁴He buffer head is operated as a 2-wire device and should ideally be driven by an AC current no greater than 1μA. A reasonable temperature estimate can be gained by driving this sensor with 10μA DC, though this is likely to cause some self-heating and also be vulnerable to thermo-electric DC offsets, particularly at higher temperatures. The five diode thermometers require excitation with currents of 10μA DC. Calibration data for all thermometer sensors are in the Excel data file that accompanies this unit.

As supplied, the heat switches require about 4V to keep them fully on at around 20 to 25 K, and will cool to the OFF state (T < 12 to 15 K) in five to ten minutes.

The cryopump heater impedances are about 200Ω for the ⁴He cryopump and about 400Ω for the ³He cryopump. During the cooling cycle it is necessary to warm the ⁴He cryopump to around 50 to 60K and the ³He cryopump to 45 to 50K. A heater current of up to 100 to 130mA or so for the ⁴He cryopump, and around 50 to 60mA for the ³He cryopump, will heat the pumps rapidly; lower heater currents will result in slower heating. Stabilisation of the pump temperatures at around 50K will typically require heater currents of around 12 to 15mA. Try to ensure that the lead-in wiring to these heaters is not unduly dissipative.

5. ATTACHING YOUR EXPERIMENT TO THE CRYOCOOLER.

This model of cryocooler provides three points at which heat may be extracted from a user's experiment mounted on a separate cold table. They are the ³He cold head, the ⁴He buffer head and the ⁴He film burner. To achieve optimum performance, only a very small load should be applied directly to the ³He cold head. The main source of cooling power is the ⁴He buffer head, which can sustain a thermal load of at least 250μW at a temperature of less than 1K. The ⁴He film burner may also be used to sink some load at around 2K.

The top surface of the ⁴He buffer head has 8 holes tapped M3 on a 40mm P.C.D. and a further axial hole tapped M4. The ³He ultra-cold head has 9 tapped holes, again M3, on a 40mm P.C.D. The film burner has 6 M3 tapped holes on the main body, in pairs on each of the three free sides. (Note: UNC threads will be substituted if requested by the customer).

While fixing experimental equipment to the cold heads, extreme care should be taken, not to torque or bend the gas pipes. Always support the cold heads against the applied torque.

The size of a GL7-type cryocooler determines its run time and temperature at a given heat loading. The cryocooler you have purchased will have been built to specified customer requirements and tested to verify that it meets its specification. Typically, the run time is limited by the buffer head, which for a standard-sized unit will last about 24 hours under the loads of 20μW on the ³He cold

head and 100 μ W on the ^4He buffer stage. Under no load, the ^3He cold head will typically run at about 260mK, the ^4He buffer head at about 865mK, and the film burner at about 1 to 1.6K. When loads are applied, the heads and film burner naturally run warmer. Load data for your specific cryocooler (which may not be typical!) are included in its accompanying Excel test data file.

5.1. Radiation shielding

The cold heads, and any cold table/experimental equipment/detector assembly you attach, must be properly radiation shielded at around 4K, in order to achieve sub-Kelvin operation. Any ancillary support structure (cold table) and experimental wiring looms must be thermally sunk to the ^4He buffer head at some point between the ^3He cold head and the 4K plate, if a satisfactory operating temperature is to be reached. Temperatures below around 300mK are only achievable if the total thermal load on the cold head is kept to a minimum (below about 20 or 25 μ W). The ^4He buffer head and film burner are designed to buffer the parasitic loads due to wiring and mechanical support structures. An additional contact to the film burner may also be desirable, though not mandatory. No other attachments to the cryocooler unit are necessary for achieve satisfactory operation.

6. OPERATION.

6.1. Pre-cool

An illustration of a typical precool to liquid nitrogen temperature in a wet dewar is shown in Figure 1. When running the cryocooler from a mechanical pre-cooler (e.g. a PT or GM cryocooler) cooling timescales will be similar unless limited by the cooling rate of the pre-cooler. The example figures in this generic manual are for a cryocooler designed to run for around 24 hours. Your own cryocooler may have been designed for a different run time and you can find similar figures, specific for your unit, in the Excel test file that accompanies your instrument.

The legend entries in the following figures are:

3-head: the ^3He cold head

4-head: the ^4He buffer head

FB: cryocooler film burner

3-Pump: the ^3He cryopump

4-Pump: the ^4He cryopump

3-Switch: the ^3He heat switch

4-Switch: the ^4He heat switch

In Figure 1 the RuO_2 sensor on the 4-head is excited with 10 μ A DC, and the effects of thermo-electric DC offsets are clear, particularly at higher temperatures. These effects are greatly mitigated at lower temperatures, and essentially vanish once the system has cooled below around 40K.





Figure 2 shows the next stage of a typical cooldown in a wet dewar during helium transfer. Timescales will be similar if using a mechanical cryocooler. The key event to watch for is the point where the heat switches turn off and the temperature of the pumps begins to rise, at around $t+25$ in the example shown. After this, active monitoring and control will be necessary to complete the cooldown (see 6.2).



6.2. Running the cryocooler

For final cooldown it is necessary to raise the temperature of both cryopumps to around 40 or 50K, and to keep their temperatures stable at that level while the heads cool to around 4 or 5K – the lower the better. In Figure 2, heater power is applied to both pumps at around t+30 and the pump temperatures are held stable at around 50K for approximately 40 minutes before proceeding to final cooldown.

A typical sequence of final cooldown events is illustrated in Figure 3. The ^4He cryopump is allowed to cool first, by turning off its heater power and turning on its heat switch. The cold head temperature and film burner temperature then start to fall rapidly. Once the ^3He cold head temperature has fallen below 2K, the ^3He cryopump is also allowed to cool by turning off its heater power and turning on its heat switch. The cold head then drops to around 250mK.

t = 5. Turn ^4He cryopump heater OFF.

t = 10. Turn ^4He cryopump heat switch ON, 4.5V.

t = 20. Turn ^3He cryopump heater OFF.

t = 30. Turn ^3He cryopump heat switch ON, 4.5V.

t = 45. Cryocooler is cooling nicely.






Figure 4 shows an example of a full cycle and run, with various applied heat loads. In this example the cryocooler ran for more than 30 hours under load. Under no applied load, the base temperatures were around 260mK for the ^3He head, 865mK for ^4He buffer head, and 1.5K for the film burner.



6.3. Typical re-cycle from cold

A re-cycle from cold is extremely simple to perform and takes less than 1 hour from start to finish, see example in Figure 5. At the start of this sequence all heater power was OFF.

- t = 3. Turn ON the ^4He cryopump heater, apply 25V (128mA, 3.2W).
- t = 5. Turn ON the ^3He cryopump heater, apply 15V (48.9mA, 733.5mW).
- t = 12. Reduce ^3He cryopump heater voltage to 4V (13.2mA, 52.8mW).
- t = 16. Reduce ^4He cryopump heater voltage to 3V (15.6mA, 46.8mW).
- t = 40. Turn ^4He cryopump heater current OFF
- t = 42. Turn ^4He cryopump switch ON, 4V.
- t = 50. Turn ^3He cryopump heater current OFF.
- t = 65. Turn ^3He cryopump heat switch ON, 4V.
- t = 80. Cryocooler is cooling nicely.



7. OPTIMISING THE PERFORMANCE OF YOUR CRYOCOOLER

It is easy to get this cryocooler to run, but it takes practice and some experimentation to achieve the best possible performance. Your particular experimental configuration will affect the thermal loadings on, and conductances between, the various parts of your cryocooler, and it may consequently alter the optimum mode of operation. You are recommended to experiment with variations on the method of operation, once some familiarity with the successful operation of the unit has been gained, in order to optimise performance for your own application.

7.1. Pre-cool to 4K.

It is most important to understand that the cold heads cool by gas convection in the gas tubes, whereas the cryopumps cool by conduction through the heat switches, while the switches are in the ON state. The heat switches will turn OFF once they cool below around 12 to 15K.

In a 'wet' cryostat, the bare unit will cool to around 80 K in 3 to 4 hours, once LN₂ is introduced into the cryostat. During L⁴He transfer, the cold heads cool rapidly while the cryopumps are warmer than about 25 or 30K, but they will cool very slowly after the cryopumps have cooled to below 25K, and the gas is adsorbed into them. If left like this, the cold heads will take two or three days to cool. The key to a rapid cooldown is to reheat both of the cryopumps to warmer than 25K and to keep them at that temperature until the cold heads have cooled.

Once the cryopumps cool to below around 25K you must wait until the corresponding heat switches have cooled sufficiently and turned OFF (i.e. to less than about 12K) before applying heater power to the cryopumps. This will warm them up to more than 25K and allow the cold heads to resume

cooling by gas convection. It will be necessary to apply heat to the cryopumps at the rate of about 50mW each, in order to keep them warm while the heads cool.

When filling a wet cryostat, if $L^4\text{He}$ is transferred at a normal rate until the main bath is full, the cryopump heat switches may not turn off until after the cryopumps are also cold. An alternate strategy is to do an initial partial transfer at a slow pace, stopping the transfer when the cryopumps reach about 20 K, and restarting an hour or two later. The optimum strategy and timings will naturally depend upon the thermal loads and masses connected to the cryocooler unit.

When running from a mechanical pre-cooler (e.g. a PT or GM cryocooler) the same considerations apply. Cooling timescales will be similar unless limited by the cooling rate of the pre-cooler.

7.1. Running the cryocooler

Once the cold heads have cooled to around 4 to 5K and both cryopump heat switches have turned OFF, the cryocooler is ready to run. The cryopumps should be warmed to around 50K to commence the cycle.

The generic method for this type of cryocooler is first to heat both pumps up to between 50 and 60 K, and to maintain them at this temperature while ensuring that the ^4He liquefaction point (in this model, the main plate) cools to below the critical point of ^4He (5.2K). In order to get the unit to operate satisfactorily, it is crucial that high liquefaction efficiency is achieved with the ^4He stage. The colder the liquefaction point gets while the pumps are warm (particularly the ^4He pump), the higher will be the liquefaction efficiency, and hence the longer the cryocooler will run before needing to be recycled. You should aim for a ^4He liquefaction temperature of below 4.6K. You should try variations on the cryopump temperatures to find a regime giving the best performance for your particular experimental set-up. Stabilising power for the ^3He cryopump should be around 40 to 50mW at 40 or 50K.

When operating from a low-powered mechanical cryocooler (e.g. 100mW @ 4K PT unit) you will probably achieve more efficient ^4He condensation by starting with the cryopumps at around 40 to 45K, rather than a higher temperature. This is because there will then be a smaller load imposed on the mechanical cryocooler due to the hot pumps. You will need to experiment in order to optimise performance for your set-up.

Once the ^4He charge is liquefied, the ^4He cryopump is allowed to cool by turning off the heater power and turning on the heat switch. This can be done more or less abruptly, depending upon the voltage applied to the switch heater. As the heat switch turns ON, the hot 4-pump imposes a large heat load onto the main plate, and this causes the main plate temperature to rise abruptly. When operating from the cold plate of a 'wet' cryostat this temperature rise is typically around 4 or 5K, to about 9K or so. When operating from a low capacity mechanical cryocooler the main plate temperature rise may briefly be greater than this, and there is the danger that if the main plate temperature should rise above about 12 or 15K the 3-pump heat switch might turn on prematurely. If this occurs then the heat switch can be turned ON slowly, by applying a lower activation voltage, and increasing it gradually, so that heat is not dumped too rapidly from the hot cryopump. If you are not in a hurry to complete the cycle you can even leave the 4-pump heat switch OFF for a while as the cryopump cools by heat leakage. It is most important to keep the 3-pump hot ($T > 35\text{K}$ or so) until the head temperatures have dropped below about 2K. The ^3He does not liquefy at the main plate, but at the film burner. The temperature of the main plate after ^4He liquefaction is not of direct

relevance to the liquefaction efficiency of the ^3He charge, although the lower you can get the head temperature before cooling the ^3He pump, the better. The precise timings and temperatures can be varied to optimise the run time for your particular application, of course.

Once the ^4He has been liquefied, and the ^4He cryopump is cooling, the cold head temperatures and film burner temperature will start to fall rapidly. The ^4He cryopump cooling curve will flatten below 25K as the charcoal begins to pump on the ^4He vapour, and the head temperature will continue to drop. You can turn off the ^3He cryopump heater power once the ^4He cryopump switch is ON, and achieve satisfactory liquefaction efficiency by waiting to turn the ^3He cryopump heat switch ON until the head temperature has dropped to below 2K or so.

Once the ^3He cryopump has begun to cool, the ^3He cold head will also cool rapidly. Final stabilisation at the operating temperature will take some time, though how long will depend upon the thermal loads that are applied to the heads by your experimental arrangement. In general, lower loads result in lower running temperatures and these require longer to achieve stabilisation. The 3-head in particular can take some while to stabilise, particularly with applied loads of less than $1\mu\text{W}$ or so. This is because the liquid ^3He has a high specific heat capacity compared to the rate at which gas evaporation at very low vapour pressure can extract latent heat. The lower the final temperature, the lower will be the corresponding saturated vapour pressure, and thus the rate at which gas evaporates.

In the example test cycle in figure 5 one can observe that the film burner temperature closely follows the ^4He cold head temperature as they cool down from 4K to around 0.9K after the 4-pump is allowed to cool. After a short while at around 0.9K the film burner temperature then starts to rise again. This signals the end of the first phase of ^4He pre-cooling. When operating the unit from a low powered pre-cooler, you are recommended to wait until this signal occurs before turning ON the ^3He cryopump heat switch. The run time of the ^4He cold head will not be greatly affected, but you may be able to achieve efficient condensation of the ^3He at a temperature below 2K without the need to raise the temperature of either cryopump as high as 50K. If the cryopump temperatures can be kept lower while still achieving efficient condensation, then the system will recycle more rapidly and also will place less strain on the pre-cooler stage. This in turn is likely to produce a more efficient and satisfactory cycle and run.

Finally, when the cryocooler is running, the heat switches for both cryopumps should be kept ON, or the cryopumps will warm up a bit and this will result in higher head temperatures.

In operation, the parasitic loading may be distributed between the ^4He buffer head and the film burner in order to optimise the ^3He head temperature or the run time. Run times of around 22 to 24 hours should be possible provided that the loads on the ^4He buffer head and the ^3He head are kept below about $150\mu\text{W}$ and $20\mu\text{W}$ respectively.

8. STANDARD PIN-OUT ASSIGNMENTS.

1	Cernox or RuO ₂	³ He HEAD	V+
14	Cernox or RuO ₂	³ He HEAD	V-
2	Cernox or RuO ₂	³ He HEAD	I+
15	Cernox or RuO ₂	³ He HEAD	I-
3	NC	NC	
4	RuO ₂	BUFFER HEAD	I+
16	RuO ₂	BUFFER HEAD	I-
5	DIODE	FILM BURNER	I+
17	DIODE	FILM BURNER	I-
6	DIODE	⁴ He PUMP	I+
18	DIODE	⁴ He PUMP	I-
7	DIODE	³ He PUMP	I+
19	DIODE	³ He PUMP	I-
8	DIODE	⁴ He PUMP SWITCH	I+
20	DIODE	⁴ He PUMP SWITCH	I-
9	DIODE	³ He PUMP SWITCH	I+
21	DIODE	³ He PUMP SWITCH	I-
10	HEATER	⁴ He PUMP SWITCH	I+
22	HEATER	⁴ He PUMP SWITCH	I-
11	HEATER	³ He PUMP SWITCH	I+
23	HEATER	³ He PUMP SWITCH	I-
12	HEATER	⁴ He PUMP	I+
24	HEATER	⁴ He PUMP	I-
13	HEATER	³ He PUMP	I+
25	HEATER	³ He PUMP	I-

Generic (i.e. standard calibration) RuO₂ sensors from Lakeshore Cryotronics are the default option on all CRC cryocoolers. Individually calibrated 'CERNOX' or RuO₂ sensors are only fitted (at additional) cost at the customer's express requirement.

Generic diode calibration curves for the cryopump and heat switch diodes are supplied as standard by CRC Ltd. A calibration curve specific to the Film Burner diode is also supplied.